



SMfit® Total Care

Installation Note

License Key:

Please use the following license key required to install SMfit® Total Care for standard display management activities:

IXXW8-VMQTB-OE6UG-NZ5KX-NQ1QP

To activate photometer support, please contact EIZO Display Technologies to request a license key update (if available for the specific photometer model).

Warning

Special care should be taken to restrict the release and distribution of SMfit® Total Care license keys and passwords. With high-level access, calibration and other device settings can be altered. Only authorized and trained users should be given this level of access. With this in mind, the license key should not be distributed and this form should be kept at an appropriately secure location.

Setup Instructions:

The following steps are necessary to install and configure SMfit® Total Care:

- Connect the displays to the workstation.
- Start the installation routine (setup.exe) and follow the setup instructions.
- Configure the workstation (e.g. enter license key, create user accounts, specify workstation properties and detect displays).
- Extract monitor data from the displays (Eizo displays only).

Please refer to the Getting Started Guide (document number A5E10000136) for detailed installation instructions.

Insert to Eizo SMfit® – 03/2008
Copyright EIZO GmbH© 2008 – All rights reserved
For additional installation or technical support:
www.eizo.eu

A5E10000241-001